

USER INSTRUCTIONS

Product code: 114.940

Available sizes: S, M, L, XL

Manufacturer: Granberg AS

Address: Bjoavegen 1442, 5584 Bjoa, Norway

A) EC Type-Examination

- I. This product is classed as Complex Design of Personal Protective Equipment (PPE) according to European PPE Directive 89/686/EEC and has been shown to comply with this Directive through the Harmonised European Standards EN420:2003 and EN374-1:2003.
- II. Notified Body responsible for certification and Article 10 compliance is BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes, MK5 8PP, UK, Notified body number 0086.
- III. Notify Body responsible for EC Quality Control System for the Final Product is BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes, MK5 8PP, UK, Notified body number 0086.



Chemical Hazards Pictogram: EN 374-1:2003; Additional information on chemical

resistance obtainable from manufacturer.

EN374:2003 permeation levels are based on breakthrough times as follows:

Performance level	1	2	3	4	5	6
Minimum Breakthrough time (mins)	10	30	60	120	240	480



Micro Organism Hazards Pictogram: EN 374-2 Penetration Test. AQL is 0,65 or EN performance level 3; Additional information obtainable from the manufacturer.

B) Performance and Limitation of Use

I. This product has been tested in accordance with EN 374-2:2003 and, EN374-3:2003 and achieved the following performance levels:

EN 374-2 Level: Passed

EN 374-3 Level:	
Sodium Hydroxide NaOH, 40%, K	6
Sulphuric Acid H ₂ SO ₄ , 96%, L	1
n-Heptane, J	2
Ethidium Bromide, 5 %	6
Formaldehyde, 37%	6
Glutaraldehyde, 50%	6
Methanol, 1.5% in water	6
Phenol, 0.1%	6

- II. Glove performance quoted is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, abrasion, degradation etc.
- III. The gloves Do Not protect against Mechanical Hazards and/or Thermal Risks.

C) Product Instruction for Use

- I. Usage – For Single Use only. If re-used, the risk of contamination and infection increases due to improper cleaning processes; and increased risk of holes and tear during re-use due to weakening of gloves by cleaning processes.
- II. Sizing – Select the right size glove for your hand.
- III. Donning – Hold glove by the bead with one hand. Align the glove thumb with your other hand thumb and slide your hand into the glove, one finger into each glove finger. Pull by the glove palm to a get a good fit. Don the other glove by the same procedure.
- IV. Inspection – Punctures or tears may occur after donning. Inspect each glove after donning, and immediately discontinue use if found damaged.
- V. During Use - Gloves polluted by blood or other organic fluids must be immediately removed as medical waste. In case of pricking or tearing, user should immediately change gloves.
- VI. In case of feeling looseness on fingers, user should change gloves.
- VII. Doffing – Hold glove bead and pull toward the finger until the glove come off.

D) Storage Condition

Store in a cool and dry place. Opened boxes should be kept away from fluorescent and sunlight.