EN: User Manual

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SV: Bruksanvisning

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Head office:
GRANBERG AS
Bjoavegen 1442, NO-5584 Bjoa, Norway
Tel +47 53 77 53 00, e-mail post@granberg.no

Swedish office / European warehouse:
GRANBERG SVERIGE AB
Schubergvägen 20, 311 74 Falkenberg, Sweden
Tel +46 (0)346 124 25, e-mail post@granberg-ab.se

granberggloves.com

ENGLISH

ART. 210.0022E – GRANBERG SURGICAL FACE MASKS TYPE IIR

Description

Surgical face masks Type IIR are made of non-woven polypropylene. Material in 3-layers with elastic loop for ears. The face masks are 100% latex free. >99% bacterial filtration efficiency.

Medical Device Regulation 2017/745. Non-sterile.
Dispose according to hospital policy and local regulation.

The face masks are made of following material:

Outer layer: Polypropylene.

Second layer: Meltblown fabric.

Inner material: Polypropylene.

Ear loops: Polyurethane.

Nose clip material: Polypropylene and steel.

Performance

This product meets the requirements of the European Standard EN 14683:2019+AC:2019: Surgical face masks, type IIR.

Test	Type I *	Type II	Type IIR
Bacterial filtration efficiency (BFE), (%)	≥ 95	≥ 98	≥ 98
Differential pressure (Pa/cm ²)	< 40	< 40	< 60
Splash resistance pressure (kPa)	Not required	Not required	≥ 16,0
Microbial cleanliness (cfu/g)	≤ 30	≤ 30	≤ 30

* Type I medical face masks should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations. Type I masks are not intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.

Face masks are in compliance with the Regulation (EC) 1935:2004 and are suitable for handling all kinds of foodstuffs for a short period of time.

Material Performance

The performance of the face masks quoted above is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, mechanical wear etc.

User guide

1. Always remember hand hygiene! Clean your hands before removing the face mask and make sure the face mask is intact.
2. Take the elastic band and place the face mask over the nose and mouth.
3. Attached the elastic band behind the ears and pull the upper edge over the back of the nose, and the lower edge under the jaw.
4. Press the bridge of the nose so that the face mask is close to the skin.
5. Make sure the edges of the face mask are close to the skin throughout the circumference.
6. When removing: The masks are removed by taking only the elastic band attached behind the ears, not the mask itself.

Replace the face mask immediately if it is difficult to breathe, it is damaged or fails to maintain a tight fit.

Storage

The face masks must be stored dry and under clean conditions where the temperature is between 10 °C - 35 °C with humidity between 45% -65%. Preferably stored in original packaging under dry conditions.

When stored as recommended, the chemical properties will remain unchanged for up to 2 years from the date of manufacture. Lifetime cannot be specified as it depends on its use. It is the responsibility of the user to determine whether the face mask is suitable for intended use.

Use limitations and Disposal

- ⚠ Components used in manufacturing may cause allergic reactions in some users. If any allergic reaction occurs, seek medical advice or visit a doctor immediately.
- EU Declaration of Conformity is accessible at: www.granberggloves.com/search
- ⚠ For single use only. Do not re-use.
- ⚠ Before usage, inspect face mask for any defect or imperfections.
- ⚠ The manufacturer cannot be held responsible for improper use of the face mask.

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